

# SPECIFICATIONS

Product Description:       **TENT**  
Part Number:               **2060**  
Style:                       **FIT TEST**

## GENERAL DESCRIPTION:

Creates a concentrated area for Qualitative Fit Testing. A saturated towelette can be clipped and hung inside the tent when using it for an Iso Amyl Acetate (*Banana Oil*) test, as well as for Bitrex or Saccharin fit tests using neubilizers.

## SPECIFICATIONS:

- Two openings for respirator adjustment.
- Hanging Strap included, to suspend tent.
- Viewing window.
- 24” Diameter supported rim for clearance around shoulder area.
- 48” High

## MATERIALS:

- Vinyl
- Nylon reinforced
- Poly coated



## WARNING

DO NOT USE FOR SMOKE FIT (Stannic Chloride) TESTING.

