

SPECIFICATIONS

Product Description: **FIT TEST KIT**
Part Number: **2040**
Style: **SACCHARIN**

GENERAL DESCRIPTION:

This kit offers qualitative OSHA compliant fit testing for disposable and reusable dust and mist respirators. Saccharin has a sweet taste to indicate failure on the respirator fit test. New disposable glass ampules make pouring solution into nebulizers easier and reduce the chance of contamination.

KIT:

- Test hood.
- Sensitivity nebulizer
- Test nebulizer
- Box of Sensitivity Solution, 6 ampules
- Box of Test Solution, 6 ampules

HOOD:

- Poly-coated
- Clear viewing window.
- Surged seams.
- Reusable

SACCHARIN:

- Glass ampule
- With pre drilled plastic cap for easy pouring.
- Plastic protection from broken glass
- One time use.
- Averages 3-5 Fit tests per tube
- Contains 2.5 cc of solution

NEBULIZER:

- Latex-free rubber air bulb for easy operation.
- Converts 5 cc (ml) of liquid into aerosol 0.3 to 5.0 microns particle size
- Complete with extra jet and replaceable plastic caps.
- Break-resistant polycarbonate plastic (Lexan®) that easily disassembles for washing
- Dimensions: 2" X 2" X 7"

AVAILABLE:

- | | | |
|---|----------------------|------------------|
| • Additional Saccharin Sensitivity solution ampules | P/N: 2040-11K | 6 per box |
| • Additional Saccharin Test solution ampules | P/N: 2040-12K | 6 per box |
| • Additional Nebulizers. | P/N: 2040-03 | Each |

